



APPLICATION

The RD-Biotech mouse immunoglobulin isotyping kit provides a rapid and easy method (one step ELISA) to characterize mouse monoclonal antibody isotypes in cell culture supernatants or purified antibodies preparations.

The kit includes ready-to-use reagents necessary to analyze 12 samples in less than 30 minutes. Buffer solutions are color coded in order to simplify pipetting steps.

PRINCIPLE OF THE ASSAY

The method employs the quantitative sandwich enzyme immunoassay technique. Anti-mouse antibodies specific to each of the common light and heavy chains are pre-coated in the wells. Samples are pipetted into microwells and Ig present in the sample are bound by the capture antibody. Then, a HRP (horseradish peroxidase) conjugated anti-mouse IgG (H+L) antibody is pipetted and incubated simultaneously with samples. After washing microwells in order to remove any non specific binding, the ready to use substrate solution (TMB) is added to microwells and color develops if the specific immunoglobulin is present in the sample. Color development is then stopped by addition of stop solution. Absorbance is measured at 450 nm.

SPECIFICITY

The method enables the detection of Light chains (Kappa and Lambda) and Heavy chains (Gamma 1, Gamma 2a, Gamma 2b, Gamma 3 and Mu) of mouse immunoglobulins.

STORAGE

All kit components are stable for 12 months when stored at 2-8°C. Do not freeze. After opening, reagents must be stored at 2-8°C, handled with care to avoid contamination and should be used within 2 months.

KIT CONTENTS

Table with 3 columns: Code, Component, Quantity. Rows include RDB3255-P (Pre-coated 8 microwells strips), RDB3255-D (Sample Diluent), RDB3255-C (Detection antibody), RDB3255-T (TMB Substrate), and RDB3255-St (Stop solution).

All the kit components are ready-to-use

ADDITIONAL MATERIAL REQUIRED

- Pipettes and tips (20-200 µl).
- ELISA plate washer (recommended)
- Microplate reader for absorbance measurements at 450 nm.
- Wash solution: H2O, 0.05% Tween 20. Other wash solutions may be used but they have to be tested with the method.

SAMPLE PREPARATION


Dilute the samples in the dilution buffer RDB3255-2 (blue solution). Cell culture supernatant: 1:20 Recommended concentration for purified Ig: 1µg/ml

Assay procedure

All steps must be performed at room temperature (RT). Bring all the reagents at room temperature 30 min before use.

Table with 2 columns: Step, Description. Steps 1-6 describe the assay procedure from adding samples to reading results.

TYPICAL RESULT



Characterization of a IgG1 antibody: yellow colour is observed in well A (corresponding to IgG1 Heavy chain), well F (corresponding to Kappa Light chain) and well H (for the positive control H+L). The well H has to be positive in order to validate the method.